

FAQs

Q: Is LIDS BY DESIGN™ safe to use?

A: Yes. The medical grade correcting strips have been eye safety tested, are latex free, and use eye area safe materials.

Q: Can I use make up after applying my LBD correcting strip?

A: Yes. For best results, the correcting strips should be applied before any makeup, oils, or moisturizers. We also recommend applying an eyelid primer over the strips before applying your makeup as usual.

Q: Is LIDS BY DESIGN™ suitable for use by men?

A: Yes. Designed for use by both men and women, the correcting strips utilize skin-like technology that is transparent and comfortable enough for all-day wear.

Q: Can I re-use the strip?

A: Each LIDS BY DESIGN™ correcting strip is intended for one time use only. Once a strip is removed, it should be discarded.

Q: How do I know what size I need?

A: In general, the heavier your eyelid, the larger the size you will need, however, eye shape may also be a factor in selecting the appropriate size. Please read our product descriptions to choose your ideal lift and fit.

Q: What if I get my placement wrong?

A: Placement takes practice. If you don't get your ideal placement on the first try, simply remove the correcting strip and re-position as needed. Once you have found the perfect placement, remove the practice strip, and replace it with a fresh strip for all day wear.

Q: How many correcting strips are in each package?

A: Each package contains 80 strips.

Q: How soon will I see a difference in my appearance?

A: Results are immediate, and will last as long as you wear the correcting strip.

Q: Does the correcting strip leave behind any residue?

A: No. When applied and removed properly, the adhesive material separates cleanly and completely from the skin.

Q: Do I have to apply make up to hide the strip?

A: No. For a nice, natural, lifted look there is no need to cover your correcting strip with make up. The strip is translucent and will blend in with most skin tones.

Q: Can I use LIDS BY DESIGN™ to correct bags under my eyes?

A: No. The correcting strips are designed and intended for use only on the upper eyelid.