

intelliWHITE™

#SmileRightNow

FUSION

2 IN 1
WHITENING SERUM

SHOPINTELLIWHITE.COM
CUSTOMER SERVICE
1-855-766-2100



intelliWHITE™

2 IN 1 WHITENING SERUM

FUSION Whitening Serum is the newest and most advanced whitening from Intelliwhite. Easy to use at home 2 in 1 whitening serum, combines 2 types of peroxide to provide both long term and short term whitening. Our patented mouth guard helps to protect your lips and gums from peroxide irritation and allows for the whitening serum to stay on teeth to effectively whiten them, no other whitening systems does this. Use 1 x a day for up to 10 days to reveal a whiter more confident smile.


SHARE YOUR RESULTS WITH US USING HASHTAG

#SmileRightNow

DIRECTIONS FOR USING FUSION 2 IN 1 WHITENING SERUM

FOR BEST RESULTS: IF YOU BRUSH YOUR TEETH PRIOR TO USING FUSION, USE A FLUORIDE-FREE TOOTHPASTE. IF YOU DO NOT HAVE A FLUORIDE FREE PASTE USE TOOTHBRUSH AND WATER

- STEP 1** Remove Whitening Guard from packaging.
Insert Whitening Guard into mouth; dry teeth well with tissue.



- STEP 2** Pump a small amount of FUSION Whitening Serum onto the whitening brush, it may take a few pumps to prime the bottle. Paint a layer of whitening serum on front surface of upper and lower teeth, keeping serum away from gum tissue.

- STEP 3** Leave on teeth for 10 Minutes

- STEP 4** Remove Whitening Guard and spit out any access. Do not eat or drink or rinse for 20 minutes, serum will continue to work. If using before bed make this your last step in your oral care routine. You can place your night appliance in your mouth after 20 minutes after whitening.

- STEP 5** Rinse Whitening Brush under water and let air dry.

Repeat above steps 1 x a day for up to 10 days.

*For optimal oral care use with IntelliWhite Smile Saver Toothpaste

*Fusion Whitening Serum can be refrigerated, we recommend refrigerating after opening to prolong the shelf life.

INTELLIWHITE ORAL CARE PRODUCTS ARE NOT A SUBSTITUTE FOR REGULAR DENTAL CARE, SO PLEASE VISIT YOUR DENTIST REGULARLY. FOR BEST RESULTS, REFRAIN FROM EATING, DRINKING OR BRUSHING FOR 20 MINUTES. INTELLIWHITE WHITENING PRODUCTS WILL ONLY WHITEN YOUR NATURAL TEETH, AND WILL ONLY WHITEN CAPS, VENEERS, CROWNS, FILLINGS AND DENTURES BACK TO THEIR ORIGINAL COLOR. DO NOT USE INTELLIWHITE-WHITENING GEL WHILE WEARING BRACES OR ANY OTHER KIND OF OVERNIGHT DENTAL APPLIANCE. DO NOT USE IF YOU ARE PREGNANT OR NURSING. PRODUCT CONTAINS XYLITOL - KEEP AWAY FROM PETS. STORE IN A COOL DRY PLACE.

INGREDIENTS: WATER, GLYCERIN, CARBAMIDE PEROXIDE, HYDROGEN PEROXIDE, TRIETHANOLAMINE, CARBOMER, XYLITOL, FLAVOR, POTASSIUM NITRATE.

DECK

TIP: During whitening treatment try to avoid drinking dark beverages or use a straw when drinking coffee, tea, or wine.